

## PRODUCT DATASHEET

Model: **PVDRNGD2**

Description: **PROTECT-A-GRID VERTICAL DR FOR FUJI D-EVO 2**

Revision: June2015

### PRODUCT ATTRIBUTES:

- Encases grid to protect against damage.
- Holds DR Panel during exams.
- Provides additional protection for DR Panel.



### MATERIAL COMPOSITION

| PART    | MATERIAL           | APPROXIMATE % |
|---------|--------------------|---------------|
| Shell   | HDP                | 98            |
| Overlay | Polycarbonate Film | 2             |

### DIMENSIONS AND PERFORMANCE

Product Dimension: 21.4 x 15.4 x 0.9 in. (54.5 x 39.1 x 2.38 cm)

Weight: 2.35 lbs. (1.07 kg)

Weight Tolerance: 400 lbs. (181.4 kg)

*Not designed for weight bearing exams. Not to be used without an image receptor.*

### FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

### BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

### STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

### REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

**Reina Imaging**

**6107 Lou Ave.**

**Crystal Lake, IL 60014 USA**

### AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148Milano Via don Gnocchi 24, Italy

+39 02 48706799

### DEVICE FAMILY CLASSIFICATION

| Country | Class                      |
|---------|----------------------------|
| USA     | FDA Class 1, 510(k) exempt |
| Canada  | Class 1                    |
| Europe  | Class 1, Annex VII, CE     |
| TGA     | Class 1                    |