

## PRODUCT DATASHEET

Model: PBHNG14

Description: **HD Mobile DR Panel Holders 2-Handle**

Revision: September 2019

**Reina**Imaging



### MATERIAL COMPOSITION

| PART    | MATERIAL               | APPROXIMATE % |
|---------|------------------------|---------------|
| Shell   | HDPE and Polycarbonate | 98            |
| Overlay | Polycarbonate Film     | 2             |

### DIMENSIONS AND PERFORMANCE

Product Dimension: 51.4 x 45 x 2.54 cm. (20.25 x 17.75 x 1.0 in)

Weight: 1 kg (2.251lb.)

Weight Tolerance: 400 lbs. (181.4 kg)

*Not designed for weight bearing exams. Not to be used without an image receptor.*

### FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

### BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

### STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

### REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

**Reina Imaging**

**6107 Lou Ave.**

**Crystal Lake, IL 60014 USA**

### AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148 Milano Via don Gnocchi 24, Italy

+39 02 48706799

### DEVICE FAMILY CLASSIFICATION

| Country | Class                      |
|---------|----------------------------|
| USA     | FDA Class 1, 510(k) exempt |
| Canada  | Class 1                    |
| Europe  | Class 1, Annex VII, CE     |
| TGA     | Class 1                    |