

PRODUCT DATASHEET

Model: CTT

Description: Clear Top Mobile Table

Revision: January 2021



Product Attributes

- **Clear polycarbonate top with scratch-resistant overlamine film**
- **4 durable swivel casters with foot locks**
- **All steel structure and welded base provides maximum strength and reliability**
- **500 lb. (227kg) patient weight limit**

Material Composition

Material	%
Steel	70%
Polycarbonate	30%

DIMENSIONS

Table Top Length: 80 in. (203.2cm)

Table Top Width: 26 in. (66.1cm)

Table Top Height: 27.9 in. (70.9cm)

Clear Span, Width: 19.25 in. (48.9cm)

Clear Span, Length: 68 in. (172.7cm)

Overall Weight: 164 lbs. (74kg)

Weight Capacity: 500 lb. (227kg)

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA