

EFFORTLESSLY PRECISE

Two-Stitch Mobile DR Panel Positioning Partner



Accomplish challenging radiographs easily in the most sensitive environments, with the Two-Stitch Mobile DR Panel Positioning Partner.

ReinaImaging


Reimagine Radiography

PRODUCT ADVANTAGES

Effortless positioning	Easy to drape, store	Easy to adjust
Supine imaging	Locking casters	Vertical lift-assist
Landscape imaging	Weighted base	Cassette compatible
Small footprint	Sturdy construction	CR or DR specific

CROSS TABLE WITH HOLDER CAUDAL AND CEPHALIC

The compact footprint of Reina's Two-Stitch Mobile DR Panel Positioning Partner makes it simple to drape and keep sterile in surgery, and easy to store. The Two-Stitch is yet another innovative and exclusive solution from Reina Imaging, your trusted provider of expertly engineered, American-made medical imaging products.



CROSS TABLE - RANGE OF MOTION

The Two-Stitch Mobile DR Panel Positioning Partner

is a multipurpose holder that allows spine exams to be performed in the most difficult imaging environments.

Two-tiered vertical travel with lift assist allows for perfect “zero-gravity” precise positioning under the surgical table, and effortless intra-operative/post-operative imaging with a single 14x17 DR panel. The Two-Stitch fits the Jackson Table clear span in the OR. Locking casters prevent unwanted movement.



CONTACT REINA IMAGING

Every customer has unique requirements. Whether you are an OEM, dealer or end user, Reina Imaging has a dedicated team that can help. Together, we'll work to understand and create solutions for your particular needs. Give us a call or send us an email today to get started with your free consultation.

ADDRESS

6107 W. Lou Ave.
Crystal Lake, IL 60014

WEBSITE

reinaimaging.com

EMAIL

customer_service@reinaimaging.com

TOLL FREE

800-752-4918

PHONE

815-356-8181

FAX

815-356-8270

